



Be Right™



QUANTAB Chloride test strips, 300 - 6000 mg/L, 0.05 - 1.0 % NaCl , 40 tests

Product #: 2751340
EUR Price: Contact Us
Available

Ideal for determining quantitative results quickly, simply dip the strip in the sample and compare the length of the colour change to a printed scale, Quantab titrators

Specifications

Measuring range:	300 - 6000 ppm Cl ⁻
Method:	Titration Strip
Number of tests:	40
Parameter:	Chloride
Platform:	Test Strip
Shipping weight:	0.2
Smallest Increments Steps:	100-200 ppm increments