



Be Right™



## Reducing agent for boiler water reagent set

**Product #:**

LCW250

Hazardous

EUR Price:

Contact Us

Available



For concentrations of 0.05-1.0 mg/L DEHA

---

### Specifications

Description:	Reducing agents for boiler water (Oxygen Scavengers)
Measuring range:	0.05 - 1.0 mg/L DEHA
Method:	Iron Reduction Method
Number of tests:	100
Parameter:	Reducing agent
Shelf life:	12 months from production date
Storage conditions:	15 - 25 °C